



Optimum® reamers

The instruction for use can be additionally requested.

Do not reuse

Before using Optimum® reamers, please see the following instruction manual.

【Contraindications and prohibitions】

Do not use in patients who are sensitive or allergic to this product.

Do not use the product unsterilized.

【Intended user】

This product shall be used only by dentists who are familiar with dental techniques.

【Target patient group】

All patients receiving dental treatment.

【Clinical benefits】

It contributes to the pulpectomy or penetration/enlargement of the root canal in dental treatment.

【Shape and materials】

- ① Working part: Stainless steel
- ② Handle: Polybutylene terephthalate(PBT)



【Intended use or effects】

Rotatable hand dental instruments used to clean the wound of the root canal or to enlarge the root canal by performing a lateral cut.

【Instructions for use】

Select instruments with an optimal size according to the intended use, and use it to debride or to enlarge and smooth the root canal wall by penetrating, cutting, and polishing with vertical reciprocating motion or lifting motion.

【Precautions for use】

1. Use the product with great care not to puncture fingers, etc. due to having a shar blade.
2. Use protective glasses, etc. to prevent eye damage. Wear a dust mask during operation to prevent sucking dust generated during operation.
3. Before use, check the product for deformation, scratch, crack, etc. outside the patient's mouth, and if any scratch or damage is found on the product, replace and discard the product.
4. If the head is slim, long, or large, avoid using the product at an unreasonable angle or under excessive pressure that may cause breaking or bending.
5. Attach a chain, thread, etc. to the opening hole of the handle to prevent fall into the oral cavity.
6. Use a rubber dam or gauze screen to prevent accidental ingestion/aspiration. For patients with a high risk of accidental ingestion/aspiration, the sitting position shall be used for treatment as much as possible, avoiding the horizontal position.
7. Long-time immersion in a corrosive solution such as EDTA and sodium hypochlorite may cause corrosion; thus, attention should be paid to such solution.

8. Do not store the product under a germicidal lamp to avoid degradation.
9. Do not use the product if the packaging is damaged.

【 Cleaning method (including treatment before sterilization)】

1. Immerse in medical detergent (enzymatic). The concentration, temperature, immersion time, etc. should be according to the instruction manual of the detergent.
2. Clean by ultrasonic with medical detergent (enzyme system). The cleaning time, etc. should be according to the instruction manual of the medical detergent.
3. When performing ultrasonic cleaning, hold it on a stand, etc. so that the working components of products do not touch each other.
4. Thoroughly rinse off the medical detergent.
5. Check visually whether there remain any adhered substances etc. If any adhered substance is found, repeat the above procedure from 1 to 4.
6. Dry the product. (Air blow, etc.)

【Sterilization method】

1. Before use, clean and sterilize by the following method.
2. Sterilization method: Put it in a sterilization tray or stand, and place it in a sterilization pack (or foil) to autoclave under either of the following conditions:

Sterilization conditions

- | | |
|---------------|---|
| Condition (1) | Temperature = 121°C, time = not less than 20 min |
| Condition (2) | Temperature = 126°C, time = not less than 15 min |
| Condition (3) | Temperature = 134°C, time = not less than 3.5 min |

3. When performing the drying process, dry at 175°C or lower for less than 30 min.

【Disposal】

1. Carefully handle when disposing of this product due to having a sharp blade.
2. This product should be discarded appropriately as medical waste.

【Distributor's name and address】

Benco Dental
295 CenterPoint Blvd., Pittston, PA 18640
Contact us : 800.GO.BENCO