



## **INSTRUCTIONS FOR USE**

### **fasTRACT® Knitted Retraction Cord (Plain)**

#### **DESCRIPTION**

- 100% cotton knitted cord designed for gingival retraction and fluid control.
- Interlocking loops enhance compressibility and placement.
- Maintains position without fraying or tearing during packing.

#### **INDICATIONS FOR USE**

- Gingival retraction prior to impression procedures
- Crown and bridge procedures
- Control of sulcular fluids
- Protection of gingival tissues during preparation

#### **INSTRUCTIONS FOR USE**

- Prepare tooth and surrounding area
- Select appropriate cord size
- Cut length to encircle tooth
- Gently pack into sulcus
- Allow 2-3 minutes
- Remove prior to procedure
- Rinse and dry area

#### **PRODUCT SIZES**

- 4006-924 – Size 00, 96 in
- 4006-933 – Size 0, 96 in
- 4006-942 – Size 1, 96 in
- 4006-951 – Size 2, 96 in



### **CONTRAINDICATIONS**

- None known when used as directed

### **WARNINGS**

- For professional dental use only
- Do not ingest
- Use caution to avoid tissue trauma

### **PRECAUTIONS**

- Single-use only; do not reuse
- Select appropriate size
- Do not force cord into sulcus
- Keep out of reach of children

### **STORAGE**

- Store in cool, dry place
- Keep container closed
- Avoid heat and moisture

RX ONLY