

BencoDental™

ValuLine® Bite Registration



- Easily stackable on the occlusal surface
- Pillow consistency allows the patient to easily close into occlusion
- Hardens into an extremely durable final set
- Tear resistant upon removal of the mouth and during trimming
- No distortion between the bite being taken and model articulation

Pack of 2 (50ml) cartridges			
Item Number	Set	Setting Time	Working Time
[5785-546]	Regular Set	2:30	1:30
[5785-555]	Fast Set	1:30	0:30
[5785-608]	Fast Set Clear	1:30	0:30
[5785-591]	Fast Set Rigid	1:30	0:30
[5785-564]	Super Fast Set	1:00	0:30

BencoDental™

For more information or to order Benco Dental products, see your Friendly Benco Rep today.

BENCO.COM
1.800.GO.BENCO