

Temporary Implant Cement



- · Dual-cured, resin based, temporary cement
- · Designed for temporary luting of prosthesis to implant abutment
- Excellent strength and durability allow for long-term cementation of implant-supported crowns
- · Radiopaque, non-eugenol, odorless and tasteless
- · Easy to blend, allowing for efficient mixing and removal
- · Automix syringe type- one syringe can be used for 100 single crowns

Kit: (1) 8g automix 1:1 cartridge, (8) mix tips, (1) Ecotip and (1) mixing pad [6290-686]

