Benco Dental

Lab Putty Shore 90



- High level of fitting accuracy and precise detail reproduction
- •90 Shore hardness
- Easy to mix with a 1:1 ratio
- •1.5 kg Base(blue) and 1.5 kg Catalyst(white)

Used for the following applications:

Fabrication matrix for hybrid over-dentures
Construction matrix for provisional restorations
Matrix for positioning anterior teeth in immediate dentures
Matrix for bar construction
Matrix for processing teeth in partial dentures
Model duplication
Preparation models for repair procedures
Insulation of teeth against plaster

Kit

[6308-489]

