



Concentrated Evacuation System Cleaner

## • Neutral pH

- · Compatible with all amalgam separators
- Dissolves and destroys a wide range of dental debris
- Fresh lemon fragrance leaves system smelling great
- Non-foaming & biodegradable
- · Easy-to-use premeasured bottle

32 oz. Bottle/Box of 100 Add 1 oz. to 32 oz. of water [5339-966]





For use in cleaning and deodorizing evacuation systems and equipment. For maximum efficiency: Clean and Deodorize with Z-Vac every other day.

 Neutral pH Non-foaming
Compatible with all amalgam separators

32 oz. [5339-966] 



For more information or to order Benco Dental" products, see your Friendly Benco Rep today.