



- Transparent medium viscosity VPS bite registration material
- Designed to be used as a matrix for fabricating standard and light cured provisional crowns and bridges
- Can be used as an occlusal template for light curing posterior composite restorations.
- Dimensionally stable, high tear strength, flexible material
- Easy to use 50 ml cartridge system

Fast Set Kit ......[4643-521]

Contains: (2) 50ml cartridges and 6 HP Teal mixing tips



To try any of our Benco Brand™ products, see your Friendly Benco Rep today. 1.800.GO.BENCO