

**BencoDental**

**BluPrint**<sup>®</sup>  
PREMIUM

## Original Formula Bite Registration



- Formulated to have a light, mousse-like consistency
- Ideal for occlusal registration, copying transfers and primary impressions
- Provides an excellent hardness/durometer with a shore A hardness of 80.
- This unique material approaches zero resistance allowing for natural jaw movement
- Non-flowing, thixotropic formula stays in place and will not drip, run or sag

**Includes 2 (50ml) cartridges, and 4 T-Mixer™ mixing tips**

Item Number	Flavor	Set	Work Time	Set Time
[5813-407]	Unflavored	Fast	30 seconds	1 minute
[5813-416]	Unflavored	Super Fast	30 seconds	30 seconds
[5813-425]	Berry	Fast	30 seconds	1 minute
[5813-443]	Berry	Super Fast	30 seconds	30 seconds
[5813-434]	Mint	Fast	30 seconds	1 minute
[5813-452]	Mint	Super Fast	30 seconds	30 seconds

**BencoDental**

For more information or to order Benco Dental products, see your Friendly Benco Rep today.

BENCO.COM  
1.800.GO.BENCO