

LED Teeth Whitening System

12% Hydrogen Peroxide



- Highly effective brush on 12% HP gel
- LED light mouthpiece promotes faster and more effective whitening and produces results in 7 to 10 days
- Includes Potassium Nitrate for sensitive teeth
- Built-in timer makes treatment convenient and precise within 15 minutes
- Each pen equals 5-10 treatments

Kit includes:

Wireless LED mouthpiece, shade guide, (3) 2ml whitening pens, owners manual, charger cord [6175-079]

Refill includes: (2) 2ml pens [6175-097]

BencoDental

For more information or to order Benco Dental products, see your Friendly Benco Rep today.

BENCO.COM 1.800.GO.BENCO